

# recommended **cleaning** & maintenance

---

Care and maintenance will ensure your dry erase surfaces are clean. Recommended products listed below.

---

## OXIVIR<sup>®</sup> TB

---



One-Step Hospital Grade Disinfectant Cleaner based on proprietary Accelerated Hydrogen Peroxide (AHP<sup>®</sup>) technology to deliver fast, effective disinfection and cleaning performance. In a simple ready to use liquid.

- Hospital Grade Disinfectant
  - Non-irritant, non-corrosive, non-bleaching, non-rinsing AHP formula
  - Virucidal, Bacteriacidal, Fungicidal, Tuberculocidal
  - Effective against 27 micro-organisms including MRSA, VRE
  - Kills viruses such as Norovirus, COVID-19, H1N1, Influenza A and Hepatitis C in 1 minute
- 

## Microban 24

---



Initially kills 99.9% of germs, including cold and flu viruses for 24 hours. To sanitize hard, non-porous, non-food contact surfaces: (1) Let stand 10 seconds, wipe clean with a damp cloth. (2) Apply and allow air to dry without wiping. Reapply if surface is scrubbed. Microban 24 does not provide 24-hour residual virus protection.

- Eliminates odors
  - Effective against *S. aureus* & *E. aerogenes* bacteria
  - Prevents mold & mildew, repeat use every 7 days
  - Effective against viruses: Influenza A H1N1, RSV and COVID-19
- 

## Isopropyl Alcohol

---



Alcohol solutions with at least 70 percent alcohol are effective against viruses on hard surfaces. First, clean the surface with water and detergent. Apply the alcohol solution (do not dilute it) and let it sit on the surface for at least 30 seconds to disinfect.

- Powerful general purpose cleaner
- Antiseptic, Antibacterial, Antiviral disinfectant and skin safe.
- Rapidly evaporates upon use
- Non-corrosive leaving zero residue